# Active<sup>3</sup><sup>®</sup>

## Erection System

Deutsch  
English  
Français  
Español  
Italiano  
Português  
Nederlands  
Dansk  
Suomi  
Norsk  
Svenska  
Polski  
Čeština  
Türkçe  
ελληνικά  
Slovenčina  
Български  
Русский  
العربية  
עברית

Medintim

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# Active<sup>3®</sup> Erection System

EN – English

Dear customer,  
Thank you for choosing our product.

Read all of the instructions for use carefully before using the product for the first time, keep them for later use and follow the instructions.

## Contents

1. General information
2. Application
- 2.1 How to prepare the application
- 2.2 How to use the erection aid
- 2.3 How to select the proper constriction ring
3. Cleaning and storing the device
4. Charging the battery
5. Side effects / contraindications
6. Disposal
7. Warranty
8. Spare and wearing parts
9. Complaints

## Delivery contents

The Active<sup>3®</sup> Erection System comprises the following parts:

- 1 pump head
- 1 transparent cylinder
- 1 large ring loader
- 1 small ring loader
- 1 cone
- 1 charging cable and a USB plug-in power supply
- 1 set of instructions
- 1 bag

## Accessories:

- 5 constriction rings in differing sizes
- 1 tube of lubricant

Please note the shelf life of the included constriction rings and the lubricant.

## 1. General information

When the penis is stiff, this is called an erection. If a man has problems getting an erection, we speak of erectile dysfunction. The Active<sup>3®</sup> Erection System is a device that helps you to get your penis stiff. Therefore, in this booklet, the Active<sup>3®</sup> Erection System is also referred to as an erection aid.

### What are the causes of erectile dysfunction?

A variety of causes can lead to erectile dysfunction. The most common ones include diabetes, thyroid disease, liver or kidney disease, high blood pressure treatment and cardiovascular disease. Prostate surgery often also leads to erectile dysfunction. Depression or exhaustion can also lead to erectile dysfunction. Sometimes, erectile dysfunction is a side effect of medication, e.g. drugs against high blood pressure. Smoking or drinking too much alcohol is another possible cause.

### How does the erection aid work?

The pump head sucks the air around the penis out of the cylinder. This creates a negative pressure (vacuum). This causes the blood vessels in the penis to fill with blood, and the penis stiffens. The constriction ring prevents the blood from flowing out. After removing the constriction ring, the blood flows back and the penis becomes soft.

## 2. Application

### 2.1 How to prepare the application

All product-related illustrations can be found on the cover flap of this booklet.

Check the device to ensure the packaging is free from external damage and complete. Make sure before use that the device and accessories are visibly undamaged and that all packaging material has been removed.

- Lay out all parts ready (Fig. 1)
- Hold the pump head as shown in Fig. 2
- Charge the pump head for at least 2 hours before using it for the first time and before a prolonged period of non-use as described in point 4

When the charging cable is connected, the motor cannot be started.

- As soon as the pump head is fully charged, you can assemble all the parts.
- Press the pump head onto the larger opening of the cylinder, as shown in Fig. 5
- First, attach the large ring loader on the other side of the cylinder and then the smaller ring loader onto it (Fig. 6)
- Insert the cone into the small ring loader (Fig. 7)
- Apply some lubricant to the tip of the cone (Fig. 8)
- First, place the largest constriction ring on the top of the cone. Slide the constriction ring over the cone onto the small ring loader. Now take off the cone (Fig. 9)
- The erection aid is now ready for use.

### 2.2 How to use the erection aid

The erection aid is for your personal use only. It may not be passed on to third parties, nor be used by them. Never use this erection aid if you have drunk alcohol or have taken recreational drugs. Do not continue using the erection aid if it does not work properly.

**Important:** Never wear the constriction ring on the penis for more than 30 minutes. Do not fall asleep with the constriction ring on the penis! Use the device only in your domestic environment, and not under the shower or in the bathtub.

- Remove the pubic hair directly on the penis and around the root of the penis.
- The device must fit the skin tightly.
- Apply some lubricant to the edge of the ring loader. (Fig. 10)
- Introduce the penis through the ring loader into the cylinder. (Fig. 11)
- Press the device tightly to your body.
- Keep the button  on the pump head permanently pressed down. (Fig. 12)
- The motor in the pump head starts and begins to suck air out of the cylinder.
- As long as you keep the button  pressed down, the device sucks air out of the cylinder.
- As soon as you release the button , the motor switches off. (Fig. 12)

### Do you see a good erection in the cylinder?

- Then push the constriction ring from the ring loader onto the stiff penis.
- Slide the constriction ring quickly and as far as possible down to the base of the penis. (Fig. 13)
- Now release the erection aid from the body. To do this, keep the button  pressed (Fig. 14). This allows air to enter the cylinder, and you can easily detach the erection aid from the body.
- The constriction ring on the penis causes the penis to remain rigid. (Fig. 15)
- Now you can have sex.

### After sex:

- Briefly pull the constriction ring apart on the flaps.
- Wait until the blood has flown back from the penis into the body.
- Then pull the constriction ring off the penis. When you use the erection aid, you will have no visible ejaculation upon orgasm (climax). (Fig. 16)
- The constriction ring retains the semen (ejaculate) in the penis.

- As soon as you remove the constriction ring, the seminal fluid flows out of the penis.
- The temporary retention of the semen in the penis is not harmful.
- After using the erection aid, clean it as described in section "How to clean the erection aid"
- Dry the parts.
- Put all parts back into the bag and keep the erection aid out of the reach of children.

### Application notes:

Application of this erection aid is easy. Nevertheless, it is better to try the erection aid initially alone and unhurried. Then you can show the erection aid to your partner. Once you both know how to use the erection aid correctly, you can enjoy sex again. Even at an advanced age, in illness or physical limitation.

### 2.3 How to select the proper constriction ring

To keep your erection as long as possible, you need to choose the proper constriction ring. If the constriction ring is too large, your penis will not remain stiff for long. If the constriction ring is too tight, it may cause pain. In case of pain, take off the constriction ring and try a larger constriction ring. When you use the erection aid for the first time, it is best to start with the largest constriction ring. If the constriction ring is too large, try a smaller one until you find the right size for yourself. The erection aid includes 2 ring loaders. The smaller one is slightly narrower, the larger one a bit wider. If the smaller ring loader is too narrow, just use the larger one. The two ring loaders are suitable for all constriction rings. Thus, it does not matter which constriction ring size you are using. If both ring loaders are too narrow, you can slip the constriction ring directly onto the cylinder. The erection aid also includes the cone. It helps you to slip the constriction ring more easily onto the ring loaders. Alternatively, you can slip the constriction ring manually onto the ring loaders without using the cone. Use only the supplied constriction rings. From time to time, the constriction rings need to be replaced. You can make repeat purchases from your dealer.

## 3. Cleaning and storing the device

Before cleaning, take the pump head off the erection aid. The pump head may only be cleaned with a damp cloth from the outside. The pump head must never be held under water as liquid can otherwise enter and damage the pump head. Do not use any detergents or solvents. All other parts can be easily cleaned with warm water and mild soap. Store the erection aid in the pouch provided in order to avoid contamination. The pump head should only be stored fully charged when not in use.

## 4. Charging the battery

Do not charge the device during use.

### How can I tell when the battery is empty?

The Active<sup>3®</sup> Erection System is equipped with a lithium ion battery. When the battery is almost empty, the button  repeatedly glows red. The device now needs to be charged. (Fig. 4 a)

- Place the charging cable on the charging contact on the pump head.
- Charge the pump head for at least 2 hours.
- During the charging process, the button  glows with a red pulsating light. (Fig. 4 b)
- When the button  glows white, the device is charged and ready for operation. (Fig. 4 c)

## 5. Side effects / contraindications

If you have a blood coagulation or blood formation disorder (e.g. sickle cell anaemia, blood cancer, bone marrow tumours, anticoagulant therapy), contact your doctor before you first use the product.

You can get a bruise (haematoma) at the base or at the top of the penis. Sometimes, there may even be slight tissue injury, when the penis was not erect for a long time. If you have diseases that cause a prolonged erection or in which sexual intercourse should be avoided, do not use the erection aid. If you have a cardiovascular disease, you should ask your doctor before you first use the product in order to clarify possible risks.

If you are not sure whether you can apply the erection aid without any problem for your health, ask your doctor or health care professional.

## 6. Disposal

### Repair and disposal information:

You must not repair the device yourself. Proper functioning will then no longer be guaranteed.

The pump head must not be opened.

In the interest of the environment, the pump head must not be disposed of with household waste at the end of its service life. Disposal may be arranged at suitable collection points in your country. Dispose of the pump head in accordance with the WEEE Guideline (Waste Electrical and Electronic Equipment). In case of questions, please contact the authority responsible for disposal.

## 7. Warranty

Almost all parts have a 2-year warranty, as prescribed by applicable law.

### The warranty does not apply:

- in case of damage resulting from incorrect use.
- for wearing parts such as constriction rings
- in the case of the customer's own fault
- once the device has been opened

## 8. Spare and wearing parts

### Use only original accessory parts.

For charging, use only the charging cable and the USB mains components supplied. You can order the replacement and wearing parts from the service address. To get the current reordering list, please contact [service@medintim.de](mailto:service@medintim.de) or your dealer.

## 9. Complaints

In case of complaints or questions regarding our service, please contact the following:

Phone: +49 6105 20 37 20

Email: [service@medintim.de](mailto:service@medintim.de)

[www.medintim.de](http://www.medintim.de)

If we ask you to send the defective product, please send the product to the following address.

KESSEL Medintim GmbH  
Nordendstr. 82–84  
64546 Mörfelden-Walldorf  
Germany

Outside Germany, please contact the dealer where you purchased the device.

Specifications			
Dimensions:	Pump head: height 87.5 mm, outer diameter 60 mm Cylinder: height 200mm, outer diameter 60 mm	Maximum negative pressure:	-350 mmHg
Weight:	900 g	Patient protection:	type BF
Power supply:	Li-Ion battery (3.7 V, 2000 mAh, 7.4 Wh)	Storage / Transport:	-25 - (+) 70°C , ≤ 90% relative humidity (non-condensing)
Protection against ingress of liquid:	IP 22	Operating environment:	(+) 5 - (+)40 °C, 15 - 90 % relative humidity
Battery lifespan:	500 charge cycles	Conformity:	IEC 60601-1: IEC 60601-1-2:2014 IEC 60601-1-11:2015 EU 93/42/EEC

### Safety and electromagnetic compatibility

When operated in accordance with the manufacturer's instructions, the Active<sup>38</sup> Erection System complies with the general requirements for the safety of medical electrical equipment in accordance with IEC Standard 60601-1 and the Supplementary Standard IEC 60601-1-2, Electromagnetic Compatibility. The battery life is 6 hours.

### Electromagnetic compatibility

The device has been tested and found to comply with the limits for medical devices pursuant to IEC 60601-1-2. These limits are designed to provide sufficient protection against harmful interference in a typical medical facility and a domestic environment.

Guidelines and manufacturer's declaration – electromagnetic immunity			
The Active <sup>38</sup> Erection System is intended for use in the following electromagnetic environment. The customer or user of the Active <sup>38</sup> Erection System should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidelines
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	Compliant	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz	Compliant	The medical electrical device is operated exclusively by battery power.
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	Not applicable (no electrical cables)	
Electrical fast transients IEC 61000-4-4	± 2 kV supply output ± 1 kV input/output power	Not applicable	Magnetic fields should be within the normal range of a typical location.
Voltage output IEC 61000-4-5	± 1 kV Differential mode ± 2 kV Common mode	Not applicable	
Mains frequency magnetic field IEC 61000-4-8	50 Hz 3 A/m (root mean square)	Compliant	
Voltage dips, short interruptions and voltage variations of the power supply input lines IEC 61000-4-11	<5% UT (> 95% UT drop) for 0.5 cycles 40% UT (60% UT drop) for 5 periods 70% UT (30% UT drop) for 25 periods <5% UT (> 95% UT drop) for 5 seconds	Not applicable	

Guidance and manufacturer's declaration – electromagnetic emissions		
The device Active <sup>38</sup> Erection System is intended for use in the following electromagnetic environment. The customer or user of the device Active <sup>38</sup> Erection System should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The medical electrical device uses RF energy only for internal functions. For this reason, the RF emissions are very low and it is not likely that they would cause interference to electronic equipment nearby.
RF emissions CISPR11	Class B	Compliant
Harmonic currents IEC 61000-3-2	Not applicable	The medical electrical device is operated exclusively by batteries.
Voltage fluctuations/ flicker emissions	Not applicable	

The Active<sup>3®</sup> Erection System meets the standards regarding electromagnetic emissions. According to current knowledge, the device is safe for use provided it is handled properly and according to the directions in these instructions for use.

Portable RF communication devices (including their accessories) should not be used at a distance of less than 30 cm from the Active<sup>3®</sup> Erection System. Failure to do so may result in reduction in the performance of the device

#### Manufacturer's Declaration – Electromagnetic Immunity II (Table 204, EN 60601-1-2)

The Active<sup>3®</sup> Erection System is intended for use in the electromagnetic environment specified below. Users of the Active<sup>3®</sup> Erection System should ensure that it is used in such an environment.

Immunity tests	IEC 60601 Test Level			Compliance Level			Electromagnetic Environment – Guidelines
	FRQ (MHz)	Modulation	Level (V/m)	FRQ (MHz)	Modulation	Level (V/m)	
Conducted RF interference according to IEC 61000-4-6 Radiated RF interference according to IEC 61000-4-3	3 V RMS 150 kHz to 80 MHz 10 V/m 80 MHz to 2.7 GHz, 1% frequency step size 80% AM @ 1 kHz			Not applicable  E1 = 10 V/m @ 80 MHz to 2.7 GHz			Interference may occur in the vicinity of equipment marked with the following symbol: 
Proximity wireless fields: IEC 61000-4-3	Test frequencies						
	FRQ (MHz)	Modulation	Level (V/m)	FRQ (MHz)	Modulation	Level (V/m)	
	385	PM 18 Hz	27	385	PM 18 Hz	27	
	450	FM	28	450	FM	28	
		5kHz devi 1kHz sine			5kHz devi 1kHz sine		
	710	PM 217 Hz	9	710	PM 217 Hz	9	
	745	PM 217 Hz	9	745	PM 217 Hz	9	
	780	PM 217 Hz	9	780	PM 217 Hz	9	
	810	PM 18 Hz	28	810	PM 18 Hz	28	
	870	PM 18 Hz	28	870	PM 18 Hz	28	
	930	PM 18 Hz	28	930	PM 18 Hz	28	
	1720	PM 217 Hz	28	1720	PM 217 Hz	28	
	1845	PM 217 Hz	28	1845	PM 217 Hz	28	
	1970	PM 217 Hz	28	1970	PM 217 Hz	28	
	2450	PM 217 Hz	28	2450	PM 217 Hz	28	
	5240	PM 217 Hz	9	5240	PM 217 Hz	9	
5500	PM 217 Hz	9	5500	PM 217 Hz	9		
5785	PM 217 Hz	9	5785	PM 217 Hz	9		

#### Manufacturer's Declaration – Recommended working clearances between the device and portable or mobile RF communications equipment (Table 206, EN 60601-1-2)

Nominal power of the emitter in W	Working clearance, depending on the transmission frequency, in m		
	150 kHz to 80 MHz $d = \{ 3.5/V1 \} \sqrt{P}$	80 MHz to 800 MHz $d = \{ 3.5/E1 \} \sqrt{P}$	800 MHz to 2.5 GHz $d = \{ 7/E1 \} \sqrt{P}$
0.01	0.04	0.04	0.08
0.1	0.11	0.11	0.22
1	0.35	0.35	0.70
10	1.11	1.11	2.22
100	3.50	3.50	7.00

For emitters, the maximum power rating of which is not specified in the table above, the recommended working clearance d in metres (m) can be determined using the equation that belongs to the respective column, where P is the maximum rated output of the emitter in watts (W) according to the manufacturer of the emitter.

The technical data is subject to change without notice for reasons of updating.

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